



*beautifully
healed skin*

BIODERMIS™
FAMILY OF QUALITY PRODUCTS

BIODERMIS Family of Scar Management Products

EPI-DERM™ Silicone Gel Sheeting

STANDARD SHEET



4.7 x 5.7 in (12 x 14.5 cm)
Small to mid-sized scars
EDGTNT-146 (1 SHEET)
EDGTNT-46 (5 SHEETS)
2 x 2.5 in. (5 x 6 cm)
Small scars
EDGTNT-123 (1 PAIR)
EDGTNT-5123 (5 PAIRS)

LARGE SHEET



11 x 16 in (28 x 40 cm)
Full-sized sheet for large scar coverage of limbs and torsos.
EDGTNT-1116 (1 SHEET)
EDGTNT-5116 (5 SHEETS)

AREOLA CIRCLES



3 x .75 in (7.6 x 1.9 cm)
Areola reconstruction procedures
EACTNT-1 (1 PAIR)
EACTNT-51 (5 PAIRS)

AREOPEXY



2.75 x 5.90 in (6.985 x 14.99 cm)
Perfect for use following the Lejour Technique. Ideal for breast reduction surgery, areola and breast reconstruction procedures.
ELPTNT-1 (1 PAIR)
ELPTNT-51 (5 PAIRS)

MASTOPEXY



12.6 x 3.8 in (32 x 9.6 cm)
Breast reduction (mammoplasty), and breast reconstruction procedures
EMPTNT-1 (1 PAIR)
EMPTNT-51 (5 PAIRS)

STRIP



1.4 x 12 in (3.5 x 30 cm)
C-section, cardiac procedures and abdominoplasty.
EDGTNT-1112 (1 PAIR)
EDGTNT-5112 (5 PAIRS)

C-STRIP



1.4 x 6 in (3.5 x 15 cm)
C-section
EDGTNT-1106 (1 STRIP)
EDGTNT-5106 (5 STRIPS)

EPITAB CIRCLES



0.75 in diameter (1.9 cm)
Preventing and managing small scar tissue.
EPTTNT-C (6 PACK)

EPITAB SQUARES



1.25 x 1.25 in (3 x 3 cm)
Squares for small scar tissue areas.
EPTTNT-S (6 PACK)

Epi-Derm™ is a soft, TNT (tacky one side only), semi-occlusive, fully encapsulating sheet developed for the prevention and management of hypertrophic and keloid scars. It is also useful for the reduction of post-burn scarring. Variety of configurations available, which may be further trimmed as needed.

EPIfoam™

EPIfoam™ offers patients and surgeons an important evolution in post-operative care. EPIfoam™ provides an excellent platform over which to apply post-operative compression garments. EPIfoam's™ unique silicone-gel surface provides a gentle adherent patient interface that helps avoid the problems associated with more aggressively adherent foam pads. EPIfoam's™ top quality, medical grade polyurethane foam helps ensure consistent compression distribution.



EPF-3 (3 PADS)
EPF-30 (30 PADS)

PRO-SIL™ Professional Silicone Therapy for Scars

Pro-Sil™ is a revolutionary, patent pending solution to scar management which provides the same industry leading technology found in Xeragel™ in a convenient, easy to apply glide-on applicator.

- All the benefits of Xeragel™
- Even easier to apply
- Compact applicator
- New formula lasts longer
- Physician tested and approved



PRO-111
PRO-511 (Box of 5)
PRO-1011 (Box of 10)

XERAGEL™ Ointment

Xeragel™ is highly recommended for use on exposed facial scars for aesthetically more appealing coverage. Xeragel™ may be worn under cosmetic foundation. Xeragel™ silicone ointment is designed to work in conjunction with Epi-Derm™ silicone gel sheeting and SilqueClenz™ low foaming cleanser as a complete scar management system.



XSO-797

EPI-NET™ Elastic Dressing

Epi-Net™ is a tubular elastic stretch net designed to serve as a secondary dressing, applying gentle pressure to keep bandages or other medical devices securely in place without adhesive tape. Ideal for securing burn or post-op dressings, hot or cold packs, IV tubes, incontinence pads and decubitus ulcer dressings. Easy to customize. Allows maximum air-flow.



Pediatric and Adult Sizes Available for:
Fingers, Toes, Wrist, Knee, Foot, Elbow, Hand, Calf, Thigh, Shoulder, Head, Chest, Abdomen, Axilla

SILQUECLENZ™/EPI-TAPE™

SilqueClenz™ is a gentle, fragrance-free, low foaming cleanser that will ensure proper care and cleaning of the gel sheeting and effective cleansing of the scar site. Epi-Tape™ is a hypo allergenic tape used to secure the gel sheeting during use.



SilqueClenz 1 fl oz (28.4 ml)
SLC-497

Epi-Tape 2 in x 10 yd (5 cm x 9.1 m)
EPT-97

SilqueClenz/Epi-Tape
SLC/EPT-97

Las Vegas, NV USA
www.biodermis.com
800-374-3376 USA
702-260-4466 International
702-260-4646 Fax
sales@biodermis.com