


## EPI-DERM NATURAL<sup>®</sup> Silicone Gel Sheeting

for keloid and hypertrophic scars

Epi-Derm<sup>®</sup> Natural is the latest silicone sheeting innovation from Biodermis. Epi-Derm Natural offers all the advantages of original Epi-Derm, but with a soft, natural colored fabric liner, designed to flex and adjust to your every move, and prevent edges from rolling under garments.

Epi-Derm<sup>®</sup>  
*Natural*  
SILICONE GEL SHEETING



C-Strip



Areola Circle



Mastopexy




Patch

### Solutions for Scars of Every Shape and Size

- **Patch** - small-sized scars
- **Standard Sheet** - mid-sized scars
- **Large Sheet** - large scars on limbs or torso
- **Large Strip** - C-section, cardiac, abdominoplasty
- **Small Strip** - small scars, knees, elbows, cannulas
- **Tabs** - circles or squares for small scars from biopsies, liposuction cannulas or mole removal
- **Areola Circles** - breast reconstruction
- **Mastopexy** - breast reduction/reconstruction
- **C-Strip** - C-section
- **Areopexy** - "lollipop" configuration used in conjunction with the Lejour Technique.

### "How does it work?"

Scars need an "ideal healing environment," meaning the appropriate balance of moisture and maximum exposure to oxygen. Epi-derm<sup>®</sup> is a fully-encapsulating sheet — it completely covers the scar treatment site. Although the entire site is covered, Epi-derm is semi-permeable, allowing oxygen to enter while keeping excess moisture out — the ideal environment for healing scar tissue.



**BIODERMIS**<sup>™</sup>  
THE SCIENCE OF SKIN  
Shop at [biodermis.com](http://biodermis.com)

1820 Whitney Mesa Drive  
Henderson, NV 89014 USA  
800-322-3729 USA  
702-260-4466 International  
702-260-4646 Fax  
[sales@biodermis.com](mailto:sales@biodermis.com)